




Product Specification Sheet

VWR Serological Pipettes, Standard Line, non-pyrogenic, sterile, 2 ml



Cat. N°	Volume (ml)	Graduations (ml)	Overfill capacity (ml)	Color code	Dimensions (mm)		
					A	B	C
612-3703	2	1/100	0,4	Green	278	2,5	6,28
612-3704	2	1/100	0,4	Green	278	2,5	6,28

Material: **100% Virgin Polystyrene**

Properties: **Non-pyrogenic**

Sterility: **Gamma Irradiated**

Accuracy: **± 3% at full volume**

Packaging 612-3703: **25 pcs/bag, 28 bags/ctn**

612-3704: **Individually wrapped, 100 pcs/bag, 8 bags/ctn**

QC Testing: **Representative production samples are collected and inspected in accordance with those product specifications**

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt, ask us.