




Product Specification Sheet

VWR Serological Pipettes, Standard Line, non-pyrogenic, sterile, 50 ml



Cat. Nr.	Volume (ml)	Graduations (ml)	Overfill capacity (ml)	Color code	Dimensions (mm)			
					A	B	C	D
612-3695	50	1	7	Black	365	5,4	8,15	18
612-3696	50	1	7	Black	365	5,4	8,15	18

Material: **100% Virgin Polystyrene**

Properties: **Non-pyrogenic**

Sterility: **Gamma Irradiated**

Accuracy: **± 3% at full volume**

Packaging 612-3695: **25 pcs/bag, 10 bags/ctn**
612-3696: **Individually wrapped, 25 pcs/bag, 4 bags/ctn**

QC Testing: **Representative production samples are collected and inspected in accordance with those product specifications**

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt, ask us.



Product Specification Sheet

VWR Serological Pipettes, Standard Line, non-pyrogenic, sterile, 50 ml
