



## Product Specification Sheet

# VWR AUTOCLAVE BAGS, BIOHAZARD

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

*VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.*

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

| Cat. No. | Description                       | Capacity (l) | L×W (mm) | Pk  |
|----------|-----------------------------------|--------------|----------|-----|
| 129-0024 | Orange with temperature indicator | 11           | 483×356  | 200 |
| 129-0025 | Orange with temperature indicator | 27           | 610×483  | 200 |
| 129-0026 | Orange with temperature indicator | 80           | 889×635  | 200 |
| 129-0027 | Orange with temperature indicator | 210          | 1143×914 | 100 |
| 129-0028 | Clear                             | 12           | 610×305  | 200 |
| 129-0029 | Clear                             | 75           | 914×610  | 200 |

**Material:** Polypropylene

**Thickness:** 50 µm

**Physical Properties:** Autoclavable up to 121°C  
High Tear resistance (ASTM D1922 tested)  
Impact Resistance (ASTM D 1709 tested)

### Other Characteristics:

All Biohazard bags have a temperature indicator patch with wording " CONTENTS HAVE BEEN AUTOCLAVED " printed with a chemically active ink. Before autoclaving the ink is clear on a white background, after autoclaving words appear in a dark brown color. Providing proof that bag has been exposed to high temperatures for sterilization.

**Disclaimer:** VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **January 2014**