



Material	103433Q
Material description	Lanthanum (III) chloride heptahydrate
Grade	AnalaR NORMAPUR
CAS Number	10025-84-0
Molecular formula	Cl ₃ La.7H ₂ O
Molecular mass	371.37

Characteristics	Specifications
Assay	Min. 98 %
Identification	Passes test
Insolubility in water	Max. 0.1 %
Ca (Calcium)	Max. 10 ppm
Cu (Copper)	Max. 1 ppm
Fe (Iron)	Max. 5 ppm
Mg (Magnesium)	Max. 1 ppm
Pb (Lead)	Max. 1 ppm
Zn (Zinc)	Max. 1 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Anja Vanhalle, Head of Laboratory - Haasrode
VWR International bvba; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.