



Material	24679.360
Material description	Paraffin oil
Grade	GPR RECTAPUR
CAS Number	8012-95-1

Characteristics	Specifications
Polycyclic aromatic hydrocarbons	Passes test
Readily carbonisable substances	Passes test
Solid paraffins	Passes test
Acidity	Max. 0.002 meq/g
Colouration	Max. 20 APHA
Density (20/4)	0.850 --> 0.885
Viscosity (40°C)	Min. 34.5 cSt
Ignition residue (SO ₄)	Max. 0.1 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250 Briare;
France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.