



Material	84530.550
Material description	n-Hexane
Grade	HiPerSolv CHROMANORM for preparative HPLC
CAS Number	110-54-3
Molecular formula	C ₆ H ₁₄
Molecular mass	86.18

Characteristics	Specifications
Assay (on anhydrous substance)	Min. 97.0 %
Acidity	Max. 0.0003 meq/g
Alkalinity	Max. 0.0002 meq/g
Evaporation residue	Max. 5 ppm
Water	Max. 50 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250 Briare;
France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.