



Certificate of Conformity

VWR DIGITAL TIMER, ONE CHANNEL WITH MEMORY AND CLOCK

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Cat.Number : 609-0136

Accuracy : $\pm 0,01$ % (3 sec.)

Colour : White

Dimensions : 2.8 in W $\times$ 2 in H
70 mm W $\times$ 51 mm H

Pk: 1

Other Characteristics:

Format of 12/24 hours

Counts down and up from 23 hours, 59 minutes, 59 seconds

Separate hour, minute, second, start/stop and reset keys

Memory function

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **July 2015**