



Certificate of Conformity

VWR DRAW ORGANISER, PS, WHITE

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

VWR Cat. No. 631-0727

Material of construction: Polystyrene

Colour: White

2 Compartment dimension: 430x80 mm

1 Compartment dimension: 430x110 mm

WxDxH : 460×307×63 mm

Other Characteristics:

- For well organised storage of various different items in drawers or at workstation. With three long compartments.
- Suitable for pipettes, thermometers, spatulas.
- Smooth surfaces with rounded corners for easy cleaning.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2016**