



Certificate of Conformity

VWR LIQUID NITROGEN CRYOGENIC FREEZER BOXES, WITH DRAIN SLOTS

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	WxDxH	Compartments	Pk
478-0247	127x125x47.6	100	1
478-0248	70x70x47.6	25	1

Material : Cardboard, White

Other Characteristics:

- Freezer boxes are constructed of durable fibreboard with protective moisture repellant coating.
- Liquid nitrogen freezer boxes feature slots for safe LN₂ drainage and reduced LN₂ consumption.
- Can accommodate dividers (sold separately) for 16 to 100 places.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **May 2016**