



# Product Specification



Material	350604T
Material description	HAEMALUM (MAYER'S) SOLUTION

Characteristics	Specifications
Suitability for microscopy(Tissue sect.)	Passes test
Nuclei	Intensive blue
Cytoplasm	Grey-blue

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International



For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.