



Certificate of Conformity

VWR MAXI BENCH PROTECTOR

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	L×W (mm)	Packed	Pk
115-0673	460x570	50 single sheets	50

Colour : White

Absorbance rate: 750 ml/m²

Other Characteristics:

- High quality, smooth, absorbent paper that quickly absorbs liquid spills with a thick, laminated polyethylene layer that prevents seeping onto the work surface.
- Retains leaked reagents in radio-chemical laboratories and prevents radioactive contamination
- Enables recovery of leaked, expensive materials
- Breakage protection for hard surfaces
- Absorption layer for water or solvent moisture chambers
- Ideal, impermeable protective material for lining animal cages

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2016**