



Certificate of Conformity

VWR SAFETY BOTTLE CARRIERS

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

VWR Cat. No.	Colour	Height W/Bail	Bucket Height	Diameter	Weight	For (ml flasks)
611-1381	Red	40 cm	27 cm	18 cm	1 kg	2500 - 4000
611-1382	Red	25.4 cm	16.5 cm	11 cm	0.20 kg	500 - 1000
611-1383	White	40 cm	27 cm	18 cm	1 kg	2500 - 4000
611-1384	Black	40 cm	27 cm	17.8 cm	1 kg	2500 - 4000
611-1385	Black	25.4 cm	16.5 cm	11 cm	0.20 kg	500 - 1000

Pack: 1

Material : Rubber

Other Characteristics:

- These high tensile strength transport carriers are designed to ensure the safe transport of bottles containing acids, alkalis and solvents.
- Handle specially shaped for carrying comfort.
- Service temperature from -45 to +110 °C - steam sterilisable.
- Locating rings in the bottom of the container ensure a gap of about 25 mm between the flask and the container wall.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2016**