



Product Specification



Material	351355C
Material description	May-Grünwald's Solution

Characteristics	Specifications
Suitability for microscopy (blood smear)	Passes test
Erythrocytes	Red
Neutrophilic granules	Faintly red
Eosinophilic granules	Deeper red
Nuclei	Blue, not violet
Lymphocyte cytoplasm	Deeper blue

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International



For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.