	Document Number	Revision	Page	Date
	TCB004	A0	1 of 1	05-09-2016
<b>PRODUCT DATA SHEET</b>		<b>Roller Bottle Cap</b>		

## 1: Material:

Cap: High Density Polyethylene Color: Purple 259


Membrane: 0.22µm Hydrophobic (Only use for Vented Cap)

## 2: Related Catalog Numbers:

VWR EU EAN. No.	VWR NA Cat. No.	Cap Style	Sterile	Packout Description
734-2821	75784-530	Vented	Yes	1/ Bag, 24/Case

## 3:Dimensions=mm

### vented cap



## 4: Features

- Vented Polyethylene cap allows consistent gas exchange but prevents passage of bacteria and fungi.
- The cap style is Screw Cap
- Vented Polyethylene cap with a 0.22µm Hydrophobic membrane
- Non-pyrogenic
- Autoclavability: No
- Shelf life: 3 years after month of production

## 5: STERILITY

- Sterilized by gamma irradiation SAL10<sup>-6</sup> (ISO11137)

## 6: PRODUCT ISO 9000 CERTIFICATION

- Manufactured in a class 100,000 room environment.
- Manufactured under EN ISO13485 :2003 and ISO9001:2008 quality management system.

Drawn	Approved
Date <b>09 May 2016</b>	Date <b>09 May 2016</b>
Name Alex Fang	Name James