



Certificate of Conformity

VWR LABORATORY ROUND BOTTLES, WITH SCREW CAP

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

VWR Cat. N°	Diameter (mm)	Height (mm)	Capacity (ml)	Thread	Pk
215-0057	186	330	5000	GL 45	1
215-0058	234	410	10000	GL 45	1
215-0059	299	505	20000	GL 45	1
215-1592	56	100	100	GL 45	10
215-1593	70	138	250	GL 45	10
215-1594	86	176	500	GL 45	10
215-1595	101	225	1000	GL 45	10
215-1596	136	260	2000	GL 45	1
215-3261	46	91	50	GL 32	10

Intended uses: Laboratory applications

Material: Borosilicate 3.3 glass

Shape: Round, with graduations, DIN GL 45 thread, pouring ring and blue PP screw cap

Compliant to: ISO 3585 Borosilicate glass 3.3 - properties
ASTM E438-92 Type 1
Eur. Pharmacopoeia 5.0, glass type I
ISO 4796 Laboratory glassware

Melamine Statement: Melamine free

Stability Data: Stable

Shelf life statement: No limit

Storage statement: Storage Conditions - common warehouse conditions.
Protected from any possible contamination.



Certificate of Conformity

Chemical Properties: Have been tested by the accredited Certification body Nr 3019 Glass Institute Hradec Kralovéd and found compliant to CSN EN ISO 3585, art.4.1,4.2,4.3a, 4 and ASTM E438-92.)
Hydrolytic class HGB 1 according to ISO 719
Hydrolytic class HGA 1 according to ISO 720

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **August 2018**