

Certificate of Conformity

VWR DISPOSABLE, ANTISTATIC WEIGHING BOATS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: **611-9174**

Lot No: **9969**

Size	I-Ø bottom mm	I-Ø top mm	Pk
Extra-Large	89	127	500

Other Characteristics:

Vacuum-molded of polystyrene homopolymer with non-toxic white pigment which is inert or minimally reactive with most common laboratory reagents. Inspected in the molding facility for cracks, incomplete or distorted shape, or other damage or defects which would compromise proper function. Sealed in polybags.

This document has been produced electronically and is valid without a signature

Date of issue: **April 2014**