



Product Specification Sheet

VWR DISPOSAL BAGS, TRANSPARENT

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat.No.	Description	Capacity(l)	L×W	Pk
129-0578	Autoclavable bags, plain	2	200×300 mm	500
129-0579	Autoclavable bags, plain	3	400×250 mm	500
129-0580	Autoclavable bags, plain	12	660×305 mm	500
129-0581	Autoclavable bags, plain	16	780×400 mm	500
129-0582	Autoclavable bags, plain	14	610×420 mm	500
129-0583	Autoclavable bags, plain	25	770×550 mm	500
129-0584	Autoclavable bags, plain	32	780×600 mm	500
129-0585	Autoclavable bags, plain	40	1100×700 mm	350
129-0586	Autoclavable bags, plain	60	1100×900 mm	100
129-0587	Autoclavable bags, with "Biohazard" label	12	660×305 mm	200
129-0588	Autoclavable bags, with "Biohazard" label	14	610×420 mm	200
129-0589	Autoclavable bags, with "Biohazard" label	32	780×600 mm	200

Material: Polypropylene

Thickness: 40 µm

Physical Properties:

- Autoclavable
- Heat resistant up to 145 °C

Appearance:

- Transparent
- No holes,dirty,scratched and air bladder



Product Specification Sheet

VWR DISPOSAL BAGS, TRANSPARENT

Anti-temperature test:

- No appearance changes and leaking after bag under 145 centigrade for 10 min

Product tensile test:

- Weld strength >12N/15mm
- Tensile strength MD/TD:25MP; Elongation:600%

Autoclave bag Artwork:

- Color code:PMS647
- Printed in 5 languages English, Deutch, French, Spanish and Italian

AQL sample plan

- Visual inspection,IL II,AQL 2.5
- Dimension inspection,IL S-3,AQL 2.5
- Function test,IL S-1,AQL 1.5

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.