



Product Specification Sheet

VWR DISPOSABLE WEIGHING AND DRYING DISHES, ALUMINIUM

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Description	Capacity (ml)	I-Ø×H (mm)	Pk
611-1350	Weighing pan	50	120×5	100
611-1351	Drying pan with hole (for solids)	30	110×11	25
611-1352	Drying Pan with Hole (for liquids)	50	110×16	25
611-1353	Weighing pan	60	100×8	50
611-1354	Weighing pan	40	70×6	50
611-1380	Weighing pan	80	70×16	100

Material: Aluminium

Intended Use: Ideal for sampling, weighing, dispensing, storage and drying applications.

Other characteristics:

- Compatible with most balances and dryers
- Oil-free
- Slightly tapered for easy stacking

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there's any doubt ask us.

Date of issue: **August 2012**

VWR International bvba/sprl, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>



Version: 2.0