



ZPACK13251000 PLASTIC BOTTLE POWDER 1 L DIN 60 UN



SPECIFICATIONS

UN-certification : 3H1/X2/S/07/D/BAM9869-M3

Material : HDPE White

Neck Type : DIN 60

Height : 97,00 mm

Length : 97,00 mm

Wide : 210,00 mm

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343