



ZPACK13252500 PLASTIC BOTTLE POWDER 2,5 L DIN 80



SPECIFICATIONS

UN-certification : 3H1/X4/S/07/D/BAM5243-M3

Material : HDPE White

Neck Type : DIN 80

Height : 280,00 mm

Wide : 127,00 mm

AVAILABLE ACCESSORIES

Key for opening : 84544.001

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343