

Material	28733.295
Material description	Trichloroethylene
Grade	GPR RECTAPUR
CAS Number	79-01-6
Molecular formula	C ₂ HCl ₃
Molecular mass	131.39

Characteristics	Specifications
Assay	Min. 97 %
IR Spectrum	Passes test
Density (20/4)	1.460 --> 1.470
Evaporation residue	Max. 50 ppm
Free alkali	Max. 0.01 meq/g
n 20/D	1.460 --> 1.480
Water	Max. 0.02 %
Cl (Chloride)	Max. 1 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250 Briare;
France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.