



Certificate of Conformity

VWR FIBREBOARD DIVIDERS

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	WxDxH	Compartments	Pk
478-0252	125x125x25.4	16	12
478-0253	125x125x25.4	25	12
478-0254	125x125x25.4	49	12
478-0255	125x125x25.4	64	12
478-0256	125x125x25.4	81	12
478-0257	125x125x25.4	100	12
478-0243	125x125x25.4	49	1
478-0244	125x125x25.4	64	1
478-0245	125x125x25.4	81	1
478-0246	125x125x25.4	100	1

Material : Cardboard, White

Other Characteristics:

- Mechanical cryogenic freezer boxes are constructed of durable fibreboard with protective moisture repellant coating.
- Liquid nitrogen (LN₂) freezer boxes feature slots for safe LN₂ drainage and reduced LN₂ consumption.
- Mechanical freezer boxes do not feature slots and cannot be used with LN₂.
- Boxes with dividers feature a lid printed with a numbered grid for superior inventory control.
- Boxes without dividers can accommodate dividers (sold separately) for 16 to 100 places.
- Fibreboard dividers can fit in any 127×127 mm LN₂ cryogenic or mechanical freezer box.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **May 2016**