



Certificate of Conformity

VWR MECHANICAL CRYOGENIC FREEZER BOXES, WITH OR WITHOUT DIVIDER

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	WxDxH	Compartments	Type	Pk
478-0239	127x127x51	-	Without dividers	1
478-0240	127x127x73	-	Without dividers	1
478-0249	127x127x47.6	81	With dividers	1

Material : Cardboard, White

Other Characteristics:

- Mechanical cryogenic freezer boxes are constructed of durable fibreboard with protective moisture repellant coating.
- Mechanical freezer boxes do not feature slots and cannot be used with LN₂.
- Boxes with dividers feature a lid printed with a numbered grid for superior inventory control.
- Boxes without dividers can accommodate dividers (sold separately) for 16 to 100 places.
- Fibreboard dividers can fit in any 127×127 mm mechanical freezer box.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **May 2016**