



Certificate of Conformity

VWR BOTTLES, SQUARE, WIDE MOUTH, WITH SCREW CAP

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	Capacity	WxDxH (mm)	Neck I-Ø	Neck Ø ext.	Pk
215-7390	100 ml	52×52×82	31 mm	40 mm	1
215-7392	250 ml	67×67×106	39 mm	49 mm	1
215-7394	500 ml	85×85×130	39 mm	49 mm	1
215-7396	1000 ml	85×85×210	39 mm	49 mm	1

Material: Bottle - High Density Polyethylene (HDPE)
Cap - Polypropylene (PP)

Other Properties:

- High resistance to non-oxidizing acids and bases.
- High chemical resistance.
- Density: 0.940 to 0.970 g/cm³
- Tamper Evident cap
- Leakproof and space saving.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **August 2016**